

RECEIVED  
CENTRAL FAX CENTER

NOV 17 2006

IN THE SPECIFICATION:

The paragraph bridging pages 7 and 8 has been amended as follows:

As PCR-amplification conditions for the ITS-DNA, the composition of the PCR reagent and the reactive conditions are shown in Table 1 and Table 2 below. The primer set may PCR amplify a sequence including ITS-DNAs or it may PCR amplify only a portion of ITS-DNAs. 1541f (~~sequence-number~~ SEQ ID NO: 1, see below) and LS23r (~~sequence-number~~ SEQ ID NO: 2, see below) are cited as examples. Other primer sets capable of PCR amplifying ITS-DNAs can also be used. After a PCR reaction, 2µl of the amplified product is confirmed using 1.5% agarose gel electrophoresis.

Please amend Table 1 as follows:

Per one reaction

Sterilized water	37.75µl
10 x PCR buffer	5.0 µl
dNTPs (2mM each)	5.0 µl (final density 0.2 mM)
<del>20 µM 1512f</del> <u>20 µM 1541f primer (sequence-number</u>	0.5 µl (final density 0.2 µM)
<u>SEQ ID NO: 1)</u>	
<del>20 µM LS23r primer (sequence-number</del> <u>SEQ ID NO: 2)</u>	0.5 µl (final density: 0.2 µM)
Taq DNA polymerase (5 units/µl) <sup>*2</sup>	0.25 µL (1.25 units)
Genome DNA	1 µl (10 ng)

\*2: HotSarTaq Polymerase by QIAGEN (Cat #203203) is used.

Page 9, line 1, has been amended as follows:

(~~Sequence-number~~ SEQ ID NO: 1)

Page 9, line 11, has been amended as follows:

(~~Sequence-number~~ SEQ ID NO: 2)